



# Straight Stairlifts






### Space Used When Seat Down



### Space Used When Seat Up



#### Measurements for Straight Stairlifts

Measurement	Bruno Elite	Brooks	Horizon Plus	Handicare 1100	Handicare Simplicity Plus Perch/1000	Brooks Sit Stand
Α	654mm	605mm	597mm to 647mm	575mm	550mm	560mm
						630mm (with grab
В	670mm	570mm	635mm	530mm	640mm	handle down)
	384mm (440mm					
	with powered					325mm (with grab
С	footplate)	310mm	340mm	280mm	370mm	handle up)
						445mm (when
D	480m to 559mm	445mm	465mm to 505mm	430mm to 520mm	675mm to 875mm	seated)

## Come and See Us!

Why not take the time to visit our showroom and try out a couple of stairlifts for yourself.

www.norfolkstairlifts.co.uk